
International Standard



1126

INTERNATIONAL ORGANIZATION FOR STANDARDIZATION • МЕЖДУНАРОДНАЯ ОРГАНИЗАЦИЯ ПО СТАНДАРТИЗАЦИИ • ORGANISATION INTERNATIONALE DE NORMALISATION

Rubber compounding ingredients — Carbon black — Determination of loss on heating

Ingrédients de mélange du caoutchouc — Noir de carbone — Détermination de la perte à la chaleur

Second edition — 1985-05-01

STANDARDSISO.COM : Click to view the full PDF of ISO 1126:1985

UDC 678.046.2 : 543.82

Ref. No. ISO 1126-1985 (E)

Descriptors : rubber industry, carbon black, tests, heating tests, determination, weight losses on heating.

Foreword

ISO (the International Organization for Standardization) is a worldwide federation of national standards bodies (ISO member bodies). The work of preparing International Standards is normally carried out through ISO technical committees. Each member body interested in a subject for which a technical committee has been established has the right to be represented on that committee. International organizations, governmental and non-governmental, in liaison with ISO, also take part in the work.

Draft International Standards adopted by the technical committees are circulated to the member bodies for approval before their acceptance as International Standards by the ISO Council. They are approved in accordance with ISO procedures requiring at least 75 % approval by the member bodies voting.

International Standard ISO 1126 was prepared by Technical Committee ISO/TC 45, *Rubber and rubber products*.

ISO 1126 was first published in 1974. This second edition cancels and replaces the first edition, of which it constitutes a minor revision.

STANDARDSISO.COM : Click to view the full PDF of ISO 1126:1985

Rubber compounding ingredients — Carbon black — Determination of loss on heating

1 Scope and field of application

This International Standard specifies a method for determining the loss on heating of carbon black for use in the rubber industry. This loss on heating is due primarily to loss of moisture, but traces of other volatile materials may also be lost.

This method is not applicable to treated carbon blacks which contain added volatile materials.

2 Principle

Heating of a test portion of carbon black for 1 h at a temperature of 105 °C or 125 °C in a weighing bottle.

Cooling of the weighing bottle in a desiccator. After weighing, calculation of the percentage loss on heating.

3 Apparatus

3.1 Oven, preferably gravity convection type, capable of maintaining a temperature of 105 ± 2 °C or 125 ± 2 °C.

3.2 Weighing bottle, squat-form, 30 mm in height and 60 mm in diameter, fitted with a ground glass stopper.

3.3 Analytical balance, accurate to ± 0,1 mg.

3.4 Desiccator.

4 Procedure

4.1 Dry the weighing bottle (3.2) and the stopper, with the stopper removed, in the oven (3.1) at a temperature of 105 ± 2 °C or 125 ± 2 °C for 30 min. Place the bottle and the stopper in the desiccator (3.4) and allow to cool to ambient temperature. Weigh the bottle with stopper to the nearest 0,1 mg.

4.2 Weigh to the nearest 0,1 mg about 2 g of carbon black into the weighing bottle.

4.3 Place the weighing bottle, test portion and stopper in the oven for 1 h at a temperature of 105 ± 2 °C or 125 ± 2 °C, with the stopper removed.

4.4 Replace the stopper and transfer the bottle and contents to the desiccator. Remove the stopper and allow to cool to ambient temperature. Replace the stopper on the weighing bottle and reweigh to the nearest 0,1 mg.

NOTE — **Precautions :**

1 Take the sample of carbon black in a tightly stoppered glass bottle or friction-top can. Allow the closed container to reach ambient temperature before starting the test.

2 Keep the weighing bottle stoppered when transferring to and from the desiccator, to prevent loss of carbon black due to air currents.

5 Expression of results

The loss on heating, expressed as a percentage by mass, is given by the formula

$$\frac{m_1 - m_2}{m_1 - m_0} \times 100$$

where

m_0 is the mass, in grams, of the weighing bottle and stopper;

m_1 is the mass, in grams, of the weighing bottle, stopper, and test portion before heating;

m_2 is the mass, in grams, of the weighing bottle, stopper, and test portion after heating.

6 Test report

The test report shall include the following particulars :

- the reference of the method used;
- the temperature used (105 °C or 125 °C);
- the results and the method of expression used;
- any unusual features noted during the determination;
- any operation not included in this International Standard, or regarded as optional.